




# CATALOG



- > Blockers
- > Animal IgG
- > Animal Serum
- > Animal Albumin
- > Secondary Antibody
- > Poly-HRP labeled Antibody



Wuhan Andigensil Biotechnology Co., Ltd. is a subsidiary of Wuhan Yuangu Biotechnology Co., Ltd., was founded in 2015. We specialize in the development and production of biological raw materials for IVD, we have self-owned animals breeding base, our core products include mouse serum, mouse IgG, bovine IgG and BSA, we can also supply other animals serum and IgG. We have got the certification of ISO13485 in 2019, and registered in TRACES of the European Commission, our products are mainly exported to the United States, Australia, Asian and European countries.



## Blockers

Product Name	Item No.	Size	Notes
YGBX01	BRX-01	1g 5g 10g	Frozen -20℃
YGBP01	BR-19	1g 5g 10g	Frozen -20℃
YGBP02	BR-20	1g 5g 10g	Frozen -20℃
YGBP03	BR-21	1g 5g 10g	Frozen -20℃
YGBM01	BR-01	1g 5g 10g	Frozen -20℃
YGBM02	BR-02	1g 5g 10g	Frozen -20℃
YGBM03	BR-03	1g 5g 10g	Frozen -20℃
YGBH02	BR-06	1g 5g 10g	Frozen -20℃
YGBG02	BR-08	1g 5g 10g	Frozen -20℃
IgG/RF Adsorbent	BR-18	1g 5g 10g	Frozen -20℃

## Animal IgG

Product Name	Item No.	Size	Notes
Mouse IgG	BMA	1g 10g 20g	Liquid, Frozen -20℃
Mouse IgG	BMA-L	1g 10g 50g	Lyophilized Powder, Storage 2-8℃
Rat IgG	BRA-L	1g 10g 50g	Lyophilized Powder, Storage 2-8℃
Guinea Pig IgG	BTA-L	1g 10g 50g	Lyophilized Powder, Storage 2-8℃

Bovine IgG	BBA-L	1g 10g 50g	Lyophilized Powder,Storage 2-8℃
Goat IgG	BGA-L	1g 10g 50g	Lyophilized Powder,Storage 2-8℃
Sheep IgG	BSA-L	1g 10g 50g	Lyophilized Powder,Storage 2-8℃
Rabbit IgG	BHA-L	1g 10g 50g	Lyophilized Powder,Storage 2-8℃
Chicken IgY	BCA-L	1g 10g 50g	Lyophilized Powder,Storage 2-8℃
BALB/c Mouse IgG	BMBA	100mg 500mg 1g	Liquid,Frozen -20℃

## Animal Albumin

Product Name	Item No.	Size	Notes
Bovine Serum Albumin (Premium Grade)	SA-001A	100g 1kg 10kg	Lyophilized Powder,Storage 2-8℃
Bovine Serum Albumin (Diagnostic Grade)	SA-001	100g 1kg 10kg	Lyophilized Powder,Storage 2-8℃
Bovine Serum Albumin (Standard Grade)	SA-001B	100g 1kg 10kg	Lyophilized Powder,Storage 2-8℃
Mouse Serum Albumin	SA-002	1g 10g 100g	Lyophilized Powder,Storage 2-8℃
Ovalbumin	SA-003	1g 10g 100g	Lyophilized Powder,Storage 2-8℃
Goat Serum Albumin	SA-005	10g 100g 1kg	Lyophilized Powder,Storage 2-8℃
Guinea Pig Serum Albumin	SA-006	1g 10g 100g	Lyophilized Powder,Storage 2-8℃
Rat Serum Albumin	SA-007	1g 10g 100g	Lyophilized Powder,Storage 2-8℃
Rabbit Serum Albumin	SA-010	10g 100g 1kg	Lyophilized Powder,Storage 2-8℃

## Secondary Antibody

Product Name	Item No.	Size	Notes
Mouse Anti-Human IgM mAb	AB-007	10mg 100mg 500mg	Frozen -20℃
Mouse Anti-Human IgG mAb	AB-030	10mg 100mg 500mg	Frozen -20℃
Rabbit Anti-Mouse IgG Antibodies	AB-036A	10mg 100mg 500mg	Frozen -20℃
Goat Anti-Rabbit IgG Antibodies	AB-037A	10mg 100mg 500mg	Frozen -20℃
Rabbit Anti-Chicken IgY Antibodies	AB-038A	10mg 100mg 500mg	Frozen -20℃
Goat Anti-Human IgG Antibodies	AB-045A	10mg 100mg 500mg	Frozen -20℃
Goat Anti-Mouse IgG Antibodies	AB-055A	10mg 100mg 500mg	Frozen -20℃
Goat Anti-Chicken IgY Antibodies	AB-056A	10mg 100mg 500mg	Frozen -20℃
Goat Anti-Human IgG Fc Antibodies	AB-057A	10mg 100mg 500mg	Frozen -20℃
Mouse anti-Human IgE mAb	AB-062A	10mg 100mg 500mg	Frozen -20℃



Rabbit Anti-Goat IgG Antibodies	AB-066A	10mg	100mg	500mg	Frozen -20℃
Goat Anti-Bovine IgG Antibodies	AB-082A	10mg	100mg	500mg	Frozen -20℃
Rabbit Anti-Rat IgG Antibodies	AB-084A	10mg	100mg	500mg	Frozen -20℃
Rabbit Anti-Guinea pig IgG Antibodies	AB-085A	10mg	100mg	500mg	Frozen -20℃
Goat Anti-Mouse IgG Antibodies	AB-055C	10mg	100mg	Cross absorbent: Human serum proteins/ Goat IgG/ Sheep IgG/ Bovine IgG/Rabbit IgG	
Goat Anti-Mouse IgG Antibodies	AB-055D	10mg	100mg	Cross absorbent: Human serum proteins/ Human IgG	
Goat Anti-Rabbit IgG Antibodies	AB-037C	10mg	100mg	Cross absorbent: Human serum proteins/ Goat IgG/ Sheep IgG/ Bovine IgG/Mouse IgG	
Goat Anti-Rabbit IgG Antibodies	AB-037D	10mg	100mg	Cross absorbent: Human serum proteins/ Human IgG	

## Animal Serum/Plasma

Product Name	Item No.	Size			Notes
Mouse Serum	SMA	1L	5L		Sterile,Frozen -20℃
Rat Serum	SRA	100ml	500ml	1L	Sterile,Frozen -20℃
Guinea Pig Serum	STA	100ml	500ml	1L	Sterile,Frozen -20℃
Rabbit Serum	SHA	100ml	500ml	1L	Sterile,Frozen -20℃
BALB/C Mouse Serum	SMBA	100ml	500ml		Sterile,Frozen -20℃
C57 Mouse Serum	SMCA	100ml	500ml		Sterile,Frozen -20℃
Mouse Plasma	PMA	500ml	1L		Sterile,Frozen -20℃
Rat Plasma	PRA	500ml	1L		Sterile,Frozen -20℃
Guinea pig Plasma	PTA	500ml	1L		Sterile,Frozen -20℃
Rabbit Plasma	PHA	500ml	1L		Sterile,Frozen -20℃

## Immunohistochemistry

Product Name	Item No.	Size			Notes
YG pen (Hydrophobic IHC pen)	PD-014	1	piece/package		/
Goat Anti-Mouse/Rabbit IgG Conjugated poly-HRP(Enhancement mode)	PD-018	100T	500T	2000T	Ready-use
Goat Anti-Mouse Conjugated poly-HRP	PD-068	100T	500T	2000T	Ready-use
Goat Anti-Rabbit IgG Conjugated poly-HRP	PD-069	100T	500T	2000T	Ready-use
Streptavidin Conjugated Poly-HRP	PD-070	100T	500T	2000T	Ready-use
Goat Anti-Mouse IgG Conjugated poly-AP	PD-071	100T	500T	2000T	Ready-use
Goat Anti-Rabbit IgG Conjugated poly-AP	PD-072	50T	200T	1000T	Ready-use



# ANDIGENSIL

Wuhan Andigensil Biotechnology Co., Ltd.



## Mouse IgG


Mouse IgG was purified by affinity chromatography, filtered by  $0.2\mu\text{m}$ , purity  $\geq 95\%$ . Mouse IgG adopts industrial production, the annual output is not less than 200kg, the batch is 5-10kg, we can meet the needs of large customers. We have industrial grade vacuum freeze-drying cabinet, our products have two dosage forms, liquid and lyophilized powder.

We can also provide other animal IgG in large quantities, including bovine, goat, sheep, rabbit, etc., all of which are not less than 95% pure.

## Mouse Serum

We have our-owned animal breeding base, all the mouse serum are got from healthy and non-immune Swiss mice, the source and quality of the serum are subject to strict quality control. Mouse serum were produced by multistage filtration process, the final filtration is  $0.45\mu\text{m}$  or  $0.2\mu\text{m}$ , the annual output exceeds 20000L. The content of IgG in the serum is not less than 9mg/ml, and the content of hemoglobin is no more than 35mg/dl.

We can also provide other animal serum, including rat, goat, sheep, rabbit, etc.



## Bovine Serum Albumin

Bovine serum albumin was purified from normal bovine serum with controlled source and quality by optimized thermal shock method and multistage filtration. The purity was  $\geq 98\%$ , monomer content  $\geq 90\%$ , and the contents of IgG, endotoxin, protease, fatty acid and other substances were strictly controlled. Products are divided into standard grade, diagnostic grade and premium grade, the annual output is more than 10 tons, we can meet the needs of customers at all levels.

We can also provide other animal albumin, including mice, rat, goat, guinea pig, chicken, rabbit, etc.



# ANDIGENSIL



[www.andigensil.com](http://www.andigensil.com)

## Contact Us

**Wuhan Andigensil Biotechnology Co.,Ltd**



No.666, Shendun 4th Road, Wuhan East Lake High Tech Development Zone, Hubei, China



+86 27 59272815

+86 15939020171 (Linda Lv)



ys\_xiong1978@163.com